



# Creating CD with Images from MediMatrix

---

*Purpose: The general steps below outline how to export DICOM images from MediMatrix to burn on a CD.*

*Requirements: The user must have access to MediMatrix as well as a CD burning program on their computer.*

## **Download the Images**

- Log into Provider Application
- Medical Records TAB> Exam Maintenance> Search for patient by name> Click  to view
- Click  to download the DICOM images to computer

## **Burn the CD**

- Find the download images on your machine> Unzip/Extract the image file
  - Find your cd burning program on your computer and drag over the images to burn
- note: no viewer is imbedded on cd, however they will open in any DICOM viewer program*

## **Need a Viewer?**

Below is a free diagnostic viewer users can download if they do not have a DICOM viewer on their computer

- MediMatrix Eye: [http://webinterstate.com/medimatrix/help/setup\\_dv.html](http://webinterstate.com/medimatrix/help/setup_dv.html)

## ***MediMatrix support:***

support@medimatrix.com

888-346-3030

www.medimatrix.com

