



Uploading EKG or Holter PDF to MediMatrix

Purpose: The instructions below explain how to upload EKG or Holter Monitor pdfs to MediMatrix allowing clients to document or send to a reading group for interpretation.

EKG and Holter Monitor Exams

- Create exam in MediMatrix – ENSURE to select correct *exam type*
- Assign technologist in orders placed screen
- Assign reading group by pushing → on exam in orders in progress screen
- Go to Daily Worklist or Medical records and click eye to open exam
- Click **UPLOAD**  in “**Upload EKG or U/S PDF Document(s) as DICOM section**”
- Find file and upload pdf (note must be acceptable file size) and SAVE
- Click **Resend Images**  in Diagnostic Images section
Exam will be routed to previously selected reading group for a report

MediMatrix support:

support@medimatrix.com

888-346-3030

www.medimatrix.com

