

# MediMatrix TM App (Android)

## **Install:**

- 1) On your Phone, navigate to the Google Play Store
- 2) Click the Search icon and type “MediMatrix”
- 3) Download the App “MediMatrix TM”
  - Supported Android version 4.4 and higher


## **Setup:**

- 1) Log into MediMatrix Provider Portal, navigate to SETUP tab> Technologists
- 2) Find your name and setup messaging as follows;
  - Use Secure Messaging = YES
  - Access Pin = setup a 5-digit access pin
  - Cancel via Secure messaging = Yes/No (if Yes selected, technician is able to cancel orders)
  - Use Tech Mobile Phone App = YES

## **Using the App:**

- 1) Dispatch will send an ‘on shift’ OR ‘exam assigned’ notification
- 2) From your phone, click the message link
- 3) User will be directed to use the App, click OPEN

(Shift Status Enabled) [https://demo.medimatrix.biz/techmobile/c?msg\\_id=l&t=Uut14YKnIPKX&z](https://demo.medimatrix.biz/techmobile/c?msg_id=l&t=Uut14YKnIPKX&z)  
On shift, please log into TechMobile.

(X-ray)  
TechMobile Access  
demo.medimatrix.biz 




# MediMatrix TM App (Android)

---

4) Upon first use, user will be prompted to choose to use location services

- Turn on Location = this will allow most accurate tracking while tracking is on (preferred)
- Close = This will disable any tracking option even when app “tracking is on”



5) Enter your 5-digit PIN, setup in previous steps, and SUBMIT

- Once logged in, users can access the App directly rather than through the text message



# MediMatrix TM App (Android)

## 6) Allowing Tracking

\*\*Due to frequency pinging limitations for devices using Version 8.0 Oreo Operating System, is preferred to run the App in the foreground for most accurate tracking\*\*

- On the BOTTOM RIGHT of App Screen user will have two options for tracking
  - TRACKING IS OFF = when displayed means tracking services are off
  - TRACKING IS ON = when displayed means tracking services are on and users can track your location
- TAP Tracking is on/Tracking is off to change tracking status

